

Microneedling Cartridges



Dr. Pen A6 ultima



Dr. Pen M8



Dr. Pen A6s ultima



Dr. Pen A11s ultima



Dr. Pen Needles

Dr. Pen needles for Microneedling therapy system (MTS)

Dr. Pen A6 ultima



Total pin	12
Application	Deep/Targeted Treatment: Acne scars, stretch marks, deeper wrinkles, hyperpigmentation.
Needle thickness	0.2mm (standard)

Dr. Pen M8



Total pin	11
Application	Targeted / Intense Treatment: Deep acne scars, wrinkles, fine lines, surgical scars, and collagen induction.
Needle thickness	0.18mm (Super Thin)



Total pin	36
Application	General Coverage / Larger Areas: Overall skin rejuvenation, hydration, mild scarring, spot removal, shrinking pores. Suitable for face and body..
Needle thickness	0.20mm (Standard)



Total pin	16
Application	Precision / General Face Treatment: Overall skin rejuvenation, fine lines, acne scarring, pigmentation, and promoting collagen.
Needle thickness	0.18mm (super thin)



Total pin	Nano
Application	Gentle Resurfacing: Enhanced absorption of serums, mild exfoliation, brightening, smoothing, and hydrating the skin. Does not puncture the dermis.
Needle thickness	Microscopic silicon pyramids



Total pin	36
Application	Broad Coverage / Larger Areas: General anti-aging, hydration, visibility of stretch marks, large pores, and treating the scalp for hair.
Needle thickness	0.20mm (standard)

Dr. Pen A6s ultima



Total pin	11
Application	Targeted / Intense Treatment: Deep acne scars, wrinkles, fine lines, surgical scars, and collagen induction.
Needle thickness	0.18mm (Super Thin)



Total pin	16
Application	Precision / General Face Treatment: Overall skin rejuvenation, fine lines, acne scarring, pigmentation, and promoting collagen.
Needle thickness	0.18mm (super thin)



Total pin	36
Application	Broad Coverage / Larger Areas: General anti-aging, hydration, visibility of stretch marks, large pores, and treating the scalp for hair.
Needle thickness	0.20mm (standard)

Dr. Pen A11s ultima



Total pin	12
Application	Targeted / Deeper Treatment: Fine lines, deep wrinkles, light scarring, acne scarring, and overall collagen induction.
Needle thickness	0.20mm standard



Total pin	24
Application	Moderate Coverage / Balanced Intensity: Versatile for various skin concerns, balancing precision and coverage.
Needle thickness	0.20mm (not specified)



Total pin	Nano round
Application	Product Absorption / Gentle Resurfacing: Enhanced absorption of serums, and general skin tone/texture improvement. Does not penetrate the dermis.
Needle thickness	Microscopic tips



DRS Derma roller

0.5mm, 1mm, 1.5mm,
2mm, 2.5mm



ZGTS Derma roller

0.5mm, 1mm, 1.5mm,
2mm, 2.5mm



Derma Stamp

