

PRP & GFC
TREATMENT PRODUCTS



Next-Level Regeneration. Real Results



PRP tubes ACD + Gel

Contains Anticoagulant citrate dextrose(ACD) and gel to separate plasma layer. Ideal for skin rejuvenation.



PRP tubes ACD + Gel + biotin


Contains Anticoagulant citrate dextrose(ACD), gel and biotin to enhance hair restoration and tissue healing process.



HA + gel PRP

Made with hyaluronic acid used fortreating osteoarthritis, especially in the knee, to relieve pain and improve function. HA provides lubrication and anti-inflammatory effects

Platelet-rich Plasma (PRP)



● Platelet-poor plasma ● Platelet-rich plasma ● Red blood cells



Golden vibrator

This is used during PRP treatment to enhance blood circulation and promote relaxation . Improves PRP efficacy results



Growth factor concentrate GFC

Professional-quality Gfc Gel Tube Designed Specifically For PRP (Platelet-Rich Plasma) procedures.



Remi centrifuge C - 852

4 x 15 ml tubes
Number of Tubes : 4 Max. Speed (RPM) : 3500
Type of Head : Swing Out Head
Supply: 220 - 240 Volts 50 Hz Single Phase

DIFFERENCE BETWEEN PRP & GFC TREATMENT

Product	GFC	PRP
Prepared from Patients Blood	Yes	No
Final Outcome	Concentrated Growth Factor derived from platelets; No RBC; No WBC	Platelet with unwanted cells like WBC and RBC
Platelet Loss	No	Yes
Complexity of process	Simple	Complex
Output Stability	7 days - at 4 C 8 hrs at room temperature	To be used immidiately
Results	Optimum and take less time	Variable and take long time
No of session required	Less	More
Pain Inflammation	Painless Treatment	Moderate Chance

