




What is GENELUX-3 ?

Principle of Q-switched Nd:YAG laser is that it generates strong energy in a very short time to instantly destroy pigment lesions in both epidermis and dermis without thermal damage to the epidermis and surrounding skin.

GENELUX-3 Advantages

- Skin Rejuvenation process will be more Gentle and Convenient, but without any consumable accessories.
- Intuitional and Simple working face, easy for operator to master the operation quickly.
- Patent special handpiece design and outlook, avoid the copy of fake companies in market.
- The newest Nano-ceramic Cavity and Insertion Extraction Device of Laser Lens, which can improve the rate of reflectivity and extend he usage life.
- The newest xenon lamp, which can extend the service life to 3 million shots.





stainless steel condenser

Adopting American stainless steel condenser,
fast heat dissipation,
constant temperature operation,
ensure the laser energy output stably.



GENELUX-3 Handpiece



1064nm-F0




Using 1064nm-F0
After applying the carbon cream, the carbon particles absorb the laser light
and explode instantly, removing impurities and keratin on the skin surface and
helping to regenerate collagen.




1064nm-F50


Before



Treatment



After




Using 1064nm-F50 for pigment disease at dermis layer with strong energy Tattoo, nevus of Ota




532nm-F50

Before/Treatment




After



Using 532nm-F50 for pigment disease at epidermis layer with strong energy Freckles, age spot, blemish



Fractional Tip (optional)




Deep skin cleaning, whitening, skin rejuvenation




Noninvasive Tattoo Removal Tip (optional)



special for eyebrow tattoo removal noninvasively



England xenon lamp, extends the service life to 3 million shots largely




Nano-ceramic cavity, improves the rate of reflectivity and extend the service life




GENELUX-3

Pigmentation Removal Principle



● superficial melanin popped out of epidermis, and dropped with the skin cuticle metabolism



GENELUX-3

Before & After



SPECIFICATION

Laser Type	Q-switched ND:YAG Laser
Mode	1064nm / 532nm
Cooling	Water cooling & Air cooling
Energy	1064nm : 50~1000mj (50mj unit) 532nm : 20~400mj (20mj unit)
Lens Type	1064-F50, 1064F-F0, 532F-50
Frequency	1~10Hz
Power Supply	220~240V, 50~60Hz
Dimension	61*45*38cm
Weight	22KG